

Allerject Instructions

Allerject is an epinephrine auto-injector, designed to be used quickly and easily in case of an allergic emergency. Allerject comes in two available doses (0.15 mg or 0.3mg) and is approximately the size of a credit card for compact handling and storing.

When the Allerject is pulled from its outer case, a built-in voice with step-by-step instructions, will automatically play to guide you through the injection process. The steps are indicated as follows:



STEP 1

Pull ALLERJECT from the outer case.

Proceed to Step 2 only when you are ready to use ALLERJECT. If you are not ready to use the device, put it back in the outer case.



STEP 2

Pull off the RED safety guard, firmly.

Do not touch the BLACK base of the auto-injector—where the needle comes out. If an accidental injection occurs, get medical help immediately.



STEP 3

Place the BLACK end against the middle of the OUTER thigh. Press firmly and hold in place for 5 seconds.

You will hear a click and hiss sound when you press ALLERJECT against your leg.

Do not inject into any other location. If you are administering ALLERJECT to a young child, hold the leg firmly in place during Step 3.

After using Allerject, you must seek medical attention. Please go to your nearest Emergency Department or call 911.

TWO AVAILABLE DOSES



ALLERJECT 0.15 mg
Intended for people who weigh 15 to 30 kg



ALLERJECT 0.3 mg
Intended for people who weigh 30 kg or more